



Flexo-Technic Handels GmbH

USED CODE 12132

Solvent recovery *Renzmann K1*



Original unit is green

Very powerfull –
1400 l Nylosolve in 24 h.

Vessel in mild steel.

For all alternative solvents.

Very easy to operate.

Saves money with the first start.

Keeps your plate-washer clean.

Renzmann K1 with vacuum
Year of make 1989 – vacuum 1996 attached very good condition - checked by Renzmann
Comes with suction pipe and flexible hose to connect with the filling pipe of the distillation vessel. Comes with emptying hose to connect with the vacuum vessel. Comes with COMPLETE documentation in English.
OPTIONAL Clean Solvent Container 500 l Used for storing clean solvent for rinsing in washing machines and for receiving the distillate from a distillation unit. * cubical container with bolted lid * all necessary connections for level switch, filling and draining * hinged lid for control of the solvent level and for sample take-out * galvanized connection pipe to the distillation unit (if the container is installed next to the distillation unit)

Technical information for the *Renzmann K1*

Unit for the recovery of high boiling solvent.



Year of manufacture	1989/1996
Volume of distillation vessel	640 litre
Constant filling volume	400 litre
Vacuum	Yes
Inline production	Possible, see chematic
Output per hour	60 - 100 litres*
Operation hours	24 hours possible
Heating up time	1 hour
Type of protection	Exe II T 3
Operation temperature	max. 200
Vacuum abs.	max. 50 mbar
Power requirement	400 V / 50 Hz / 3Ph
Power consumption oil-heater	22,5 kW
Consumption cooling water	0,8 m ³ /h
Payment	50% with your order 50% before crating.
Price, ex-works, un crated	EURO 22.000

* Depends on type of solvent, type and percentage of contamination

Technical changes are subject to alteration.